



**Omega  
Fatty Acid  
Isomers**

**Separation of Native and Unsaturated  
Fatty Acids**


**Author: AB  
Date:  
i-Mail-No.**

Omega-3- and Omega-6- fatty acids have an important impact on health with positive effects on heart diseases and arteriosclerosis being reported. These fatty acids are essential, as they cannot be synthesised by humans. For this reason an adequate supply via nutrition is mandatory. Often they are applied as dietary supplements or added to conventional food products. An appropriate ratio of Omega-3 and Omega-6 is recommended.

Therefore, a reproducible and simple analysis of fatty acids in the field of food industry is fundamental. The shown application enables the successful separation of structurally very similar compounds – **without derivatisation of fatty acids!**



Well known omega fatty acids are  $\alpha$ -linolenic acid (C18:3  $\alpha$ , Omega-3) and  $\gamma$ -linolenic acid (C18:3, Omega-6). They differ only in their position of the double bonds.






**Omega  
Fatty Acid  
Isomers**

**Separation of Native and Unsaturated  
Fatty Acids**

**Author: AB  
Date:  
i-Mail-No.**



<b>Column</b>	YMC-Triart C18 ExRS TAR08S03-1546PTH 3 µm particle size 8 nm pore size 150 x 4.6 mm ID	
<b>Mobile Phase</b>	A: H <sub>2</sub> O + 0.5 % H <sub>3</sub> CCOOH B: ACN + 0.5 % H <sub>3</sub> CCOOH	
<b>Gradient</b>	Time [min]	% B
	0	76.5
	22	78.5
	27	78.5
	27.1	90
	35	90
<b>Flow rate</b>	1 mL/min	
<b>Detection</b>	205 nm	
<b>Temperature</b>	30 °C	
<b>Substances</b>	1: Eicosapentaenoic acid (C20:5; 0.5 mM) 2: Alpha-Linolenic acid (C18:3 α; 0.5 mM) 3: Gamma-Linolenic acid (C18:3 γ; 0.5 mM) 4: Docosahexaenoic acid (C22:6; 0.5 mM) 5: Palmitoleic acid (C16:1; 3 mM) 6: Arachidonic acid (C20:4; 0.5 mM) 7: Linoleic acid (C18:2; 0.5 mM) 8: Oleic acid (C18:1; 3 mM)	
<b>Injection</b>	20 µL	



**Omega  
Fatty Acid  
Isomers**

**Separation of Native and Unsaturated  
Fatty Acids**

**Author: AB  
Date:  
i-Mail-No.**

With YMC-Triart C18 ExRS you achieve a separation of fatty acids:

- Scalable from U(HPLC) to HPLC
- Without derivatisation of the analytes
- Suitable for structurally similar compounds

**Ask for your YMC-Triart C18 ExRS free trial column now!**

For further information please contact:

Dr. Anja Bräutigam  
phone: +49 (0)2064 427-283  
E-mail: [braeutigam@ymc.de](mailto:braeutigam@ymc.de)

