





UHPLC Size Exclusion Columns

YMC-Pack Diol 2 μm




YMC-Pack Diol UHPLC

- Faster analysis and increased throughput
- Excellent reproducibility
- Improved resolution
- Long-term stability
- Easy method transfer from HPLC to UHPLC



	YMC-Pack Diol-200	YMC-Pack Diol-300
Base particle	Silica	
Functional group	Dihydroxypropyl	
Particle size	2 µm	
Pore size	20 nm / 200 Å	30 nm / 300 Å
Recommended pH range	5.0-7.5	
Pressure limit	45 MPa (6,500 psi)	
Molecular weight range	5,000 – 300,000 Da	20,000 – 1,000,000 Da


Higher resolution for analysis of monoclonal antibodies



All three particle sizes show identical separation patterns for monoclonal antibody analysis. This allows easy method transfer between HPLC and UHPLC. A method developed using conventional HPLC can be directly transferred to UHPLC using a 2 µm YMC-Pack Diol column.

YMC-Pack Diol UHPLC columns greatly improve the resolution between aggregates and the monomer peak. Furthermore, a shoulder peak which elutes after the monomer peak can partially be separated using the 2 µm column.


High throughput analysis of monoclonal antibodies



By using a 150 mm length column, 50% shorter run time can be achieved with the same resolution as for a 5 µm 300 mm length column. Doubling the throughput can be achieved.

The backpressure was only 6.6 MPa even for the 300 mm column. Therefore, YMC-Pack Diol 2 µm columns can be used with both UHPLC and HPLC.


Excellent batch to batch reproducibility



YMC-Pack Diol UHPLC columns have excellent batch-to-batch reproducibility. This makes YMC-Pack Diol 2 µm columns the ideal choice for the quality control of bio-based drugs including monoclonal antibodies.

YMC-Pack Diol UHPLC

Long-term stability



YMC-Pack Diol UHPLC columns maintain their performance for more than 500 injections of sample during monoclonal antibody analysis. This ensures reproducible and reliable quality control of bio-based drugs including monoclonal antibodies.

Ordering information

Phase	Particle size [μ m]	Pore size [nm]	Column ID [mm]	Column dimension Length [mm]	Part number
YMC-Pack Diol-200	2	20	4.6	150	DL20S02-1546PTH
				300	DL20S02-3046PTH
YMC-Pack Diol-300	2	30	4.6	150	DL30S02-1546PTH
				300	DL30S02-3046PTH

YMC Co., Ltd.

YMC Karasuma-Gojo Bld. 284 Daigo-cho,
Karasuma Nisuru Gojo-dori Shimogyo-ku,
Kyoto 600-8106 Japan
TEL. +81(0)75-342-4515, FAX +81(0)75-342-4550
www.ymc.co.jp

YMC Europe GmbH

Schöttmannshof 19
D-46539 Dinslaken
Germany
TEL. +49(0)2064/427-0, FAX +49(0)2064/427-222
www.ymc.de

YMC America, Inc.

941 Marcon Boulevard Suite 201
Allentown, PA18109 USA
TEL. +1-610-266-8650, FAX +1-610-266-8652
www.ymcamerica.com

YMC Taiwan Co., Ltd.

1F., No. 76, Ln. 191, Sec. 3, Minquan E. Rd.,
Songshan Dist., Taipei City 105,
Taiwan (R.O.C.)
TEL. +886-3-2150-630, FAX +886-3-2150-286
www.ymctaiwan.com

YMC Korea Co., Ltd.

310 Woolim W-City, 618 Sampyeong-dong,
Bundang-gu, Sungnam-si, Gyeonggi-do
463-400 Korea
TEL. +82-31-716-1631, FAX +82-31-716-1630
www.ymckorea.com

YMC India Pvt. Ltd.

A-154, Eros Boulevard Hilton Hotel Complex Plot
No.13-B, District Centre Mayur Vihar Phase-I
New Delhi-110091 India
TEL. +91-11-45041601;45041701, FAX +91-11-45041901
www.ymcindia.com

YMC Co., Ltd. Shanghai Rep. Office

Far East International Plaza A2404
No. 319 Xianxia Road, Shanghai 200051
P.R. China
TEL. +86-21-6235-1388, FAX +86-21-6235-1398
www.ymchina.com

YMC Singapore Tradelinks Pte. Ltd.

6 Indus Road, #13-09 Emerald Park, Tower 1,
Singapore 169588
TEL. +65-9023-7617
www.ymc.sg