

## Purification of Peptides with Full Flexibility

Purification is the most critical step in the manufacturing process of **peptide therapeutics**. The right choice of chromatography media is crucial for **cost-effective production**. With its **wide pH range** (pH 2-10), YMC-Triart Prep C18-S provides you with **full flexibility** in the method development of peptide purification. **Simple scale-up** procedures ensure the reproducible result at manufacture-scale. A method for the purification of liraglutide with high resolution (antidiabetic peptide therapeutic, marketed by Novo Nordisk as Victoza®.) was successfully developed with YMC-Triart Prep C18-S under alkaline condition. The purity obtained for the target compound was 99.5%.




Figure 1:  
The purification of liraglutide with YMC-Triart Prep C18-S (10  $\mu$ m, 12 nm, 250 x 10 mm ID).

Column: YMC-Triart Prep C18-S (10  $\mu$ m, 12 nm, 250 x 10 mm ID)  
 Eluent: A) 20 mM HCOONH<sub>4</sub>-NH<sub>3</sub> (pH 8.5)  
           B) Acetonitrile  
 Gradient: 30% - 50% B (0 - 50 min)  
 Flow rate: 4.7 mL/min  
 Temperature: Ambient  
 Detection: UV at 215 nm  
 Injection: 3 mL (Crude 20.0 mg/mL)

## Purification of Peptides with Full Flexibility

### Improved Resolution under alkaline condition

During the scouting process of liraglutide purification, under alkaline condition (pH 8.5), the chromatogram shows an obvious better separation of the peaks with a new impurity peak also appearing.




Figure 2:  
Optimization of pH for the purification of liraglutide.

Column: YMC-Triart C18 (3  $\mu$ m, 12 nm, 100 x 3.0 mm ID)  
 Eluent: A) Different buffer showed in the figure B) Acetonitrile  
 Gradient: 40% - 75% B (0 - 10 min)  
 Flow rate: 0.43 mL/min  
 Temperature: 35°C  
 Detection: UV at 215 nm  
 Injection: 6  $\mu$ L (Crude 0.5 mg/mL) = 3  $\mu$ g loading

## Purification of Peptides with Full Flexibility

### Easy Scale-up with YMC-Triart Prep

The developed method with YMC-Triart Prep can easily be scaled-up. Below is an example of the theoretical scale-up calculation for the developed method for liraglutide purification with YMC-Triart Prep C18-S. With an YMC-Triart C18-S (250 x 600 mm I.D.) column, up to 800 g liraglutide can be purified per day.

Table 1: Scale-up calculations for liraglutide purification.

<b>Column</b>	YMC-Triart C18-S (10 µm, 12 nm)		
<b>Eluent</b>	A) 20 mM HCOONH <sub>4</sub> -NH <sub>3</sub> (pH 8.5)		B) Acetonitrile
<b>Gradient</b>	30-50% B (0-50 min)		
<b>Detection</b>	UV at 215 nm		
<b>Temperature</b>	Ambient		
<b>Cycle time</b>	60 min/run - 8 cycles/day		
Column dimension	250 x 100 mm ID	250 x 450 mm ID	250 x 600 mm ID
Flow rate	0.47 L/min	9.52 L/min	16.92 L/min
Loading / run	6.0 g	121.5 g	216.0 g
Fraction volume /run	1.4 L	28.6 L	50.8 L
Liraglutide recovery / run	2.6 g	53.4 g	94.9 g
Liraglutide recovery / day	20.8 g	427.2 g	759.2 g

## Conclusions

### Benefits of YMC-Triart Prep for liraglutide purification:

- An optimized method at high pH with improved resolution
- Up to 4-fold longer lifetime than conventional silica materials
- High loadability and high productivity
- Easy scale-up procedures

