

# Reversed Phase Method Development



Step 1

## Characteristics

### Applicational characteristics

- HPLC/UHPLC-method?
- Which detector? LC/MS?
- Isocratic/gradient?

### Analyte characteristics

- Hydrophobicity, polarity, ionicity
- Structure/molecular weight
- Stability
- How do the analytes differ?
- Matrix

Step 2

## Screening

### Typical screening conditions:

- Steep gradient
- Mostly short column e.g. 50 mm
- ID depends on compound, pressure

### Columns

#### Screening kit

1. YMC-Triart C18
2. YMC-Triart C18 ExRS
3. YMC-Triart C8
4. YMC-Triart Phenyl
5. YMC-Triart PFP

X

### Solvent

- 1 Acetonitrile
- 2 Methanol


Gradient  
5–90% solvent

X

### pH (solution)

- I. Acidic
- II. Neutral
- III. Basic

### Choose the most promising conditions



Example 6 pigments on YMC-Triart columns  
Column: 50 x 2.0 mm ID  
Gradient: 5–90%B(0–5 min), 90%B(5–7 min), 5%B(7–12 min)  
Flow rate: 0.2 mL/min  
Temperature: 40 °C  
Detection: UV at 254 nm

1. Acid green 16
2. Acid blue 1
3. Acid red 52
4. Acid blue 3
5. Methyl orange
6. Methyl red

Step 3

## Optimisation

### Adjustment options


- Mobile phase ratio
- Gradient slope/isocratic
- Precision tuning of temperature
- Column dimensions
- Particle size
- Add org. modifier
- pH (additives)

5–90%B(0–20 min)

Separated by decreasing in gradient slope

0–15%B(0–8 min)  
15–40%B(8–15 min)  
40–90%B(15–20 min)

### Optimisation of mobile phase ratio and gradient slope



Example cold medicine on Hydrosphere C18  
Column: Hydrosphere C18 (5 µm, 120 Å) 150 x 4.6 mm ID  
Eluent: A) 20 mM phosphate buffer (pH 2.5)  
B) methanol  
Flow rate: 1.0 mL/min  
Temperature: 37 °C  
Detection: UV at 210 nm (0–15 min), 235 nm (15–25 min)

1. Thiamine hydrochloride
2. Unknown
3. L-Ascorbic acid
4. Maleic acid
5. di-Methylephedrine hydrochloride
6. Dihydrocodaine phosphate
7. Saccharin sodium
8. Caffeine
9. Chlorpheniramine
10. Ibuprofen

### Buffer Selection

Buffers should be chosen as following: **pH 2 values under pKa** of the analyte for acidic compounds and for **basic** compounds **2 values above the pKa**. For LC/MS methods it is recommended to use buffer concentrations <15 mM.

Buffer	pKa	Buffer range [pH]	Standard Concentration	LC/MS compatibility
Trifluoroacetic acid (TFA)	<1.0	–	0.01–0.1%	✓
Phosphoric acid	2.1	–	0.01–0.1%	✗
Ammonium dihydrogen phosphate (Na <sup>+</sup> , K <sup>+</sup> salt)	2.1	1.1–3.1	5–50 mM (<20 mM recommended)	✗
Formic acid	3.7	–	0.1–1.0%	✓
Ammonium formate (Na <sup>+</sup> , K <sup>+</sup> salt)*	3.7	2.7–4.7	5–50 mM	NH <sub>4</sub> <sup>+</sup> salt: ✓ [Na <sup>+</sup> , K <sup>+</sup> salt: ✗]
Acetic acid	4.8	–	0.5–5.0%	✓
Ammonium acetate (Na <sup>+</sup> , K <sup>+</sup> salt)*	4.8	3.8–5.8	5–50 mM	NH <sub>4</sub> <sup>+</sup> salt: ✓ [Na <sup>+</sup> , K <sup>+</sup> salt: ✗]
Ammonium hydrogen phosphate (Na <sup>+</sup> , K <sup>+</sup> salt)	7.2	6.2–8.2	5–50 mM (<20 mM recommended)	✗
Triethylamine acetic acid (TEAA)	–	4.6–6, 10–11	<20 mM	✓
Ammonium formate, ammonium acetate**	9.2	8.2–10.2	<20 mM	✓
Sodium phosphate, potassium phosphate	12.3	11.3–13.3	<10 mM	✗
Ammonium bicarbonate**	–	8.5–10.5	<10 mM	✓

\*adjusted with acid \*\*adjusted with ammonia

### YMC-Triart Phase Specifications

	C18	C18 ExRS	Bio C18	C8	Bio C4	Phenyl	PFP
Base	organic/inorganic silica						
Stationary phase	C18 (USP L1)	C18 (USP L1)	C18 (USP L1)	C8 (USP L7)	C4 (USP L26)	Phenyl (USP L11)	Penta-fluorophenyl (USP L43)
Particle size	1.9, 3 and 5 µm						
Pore size	12 nm	8 nm	30 nm	12 nm	30 nm	12 nm	12 nm
Carbon content	20%	25%	—	17%	—	17%	15%
Endcapping	multi-stage	multi-stage	multi-stage	multi-stage	multi-stage	multi-stage	none
pH range	1–12	1–12	1–12	1–12	1–10	1–10	1–8
Temperature range	pH<7: 90 °C pH>7: 50 °C	pH<7: 90 °C pH>7: 50 °C	pH<7: 90 °C pH>7: 50 °C	pH<7: 90 °C pH>7: 50 °C	pH<7: 90 °C pH>7: 50 °C	50 °C	50 °C
100% aqueous eluents	✓	✗	✓	✗	✓	✓	✓