



Clindamycin hydrochloride

Column : YMC-Triart C18 (5 μ m, 12 nm)
 250 X 4.6 mm I.D.
 Eluent : 50 mM KH_2PO_4 (pH 7.5 adjusted by 8 M KOH)/acetonitrile (55/45)
 Flow rate : 1.1 mL/min
 Temperature : 25°C
 Detection : UV at 210 nm
 Injection : 20 μ L (1.0 mg/mL)

(第十五改正日本薬局方第二追補記載条件；定量法)