



Glass Column Packing Service

YMC
EUROPE GMBH



Save time with prepacked YMC glass columns!

Are you lacking the time to pack your column?
Or are you looking for a ready to use product?

Here is the solution!

Rely on the experience of YMC and ask us for prepacked glass columns. Our packing expertise is your advantage.

You can choose from the following services:

Packing of glass columns for your IEX, AC, SEC, HIC and RP application

- Inner diameter 5–80 mm
- Length up to 1,000 mm
- Customised dimensions possible.

Packing of new columns or refilling of your existing columns

Of course you can order new prepacked columns.

Or do you need a refill of your existing YMC columns? The refill can be combined with additional maintenance services, e.g. changing of frits and O-rings.

Packing of stationary phases from YMC or from other stationary suppliers

Is your method restricted to the use of a specific stationary phase? Of course, we can pack your YMC glass column according to your application with the phase of your choice.

If the stationary phase is not manufactured by YMC, please provide it for packing service.

How it works:



Please contact us for further information and inquiries:

Dr. Mathias Hehn
Head of Laboratories
Phone: +49 (0)2064 427-283
Fax: +49 (0)2064 427-222
Email: hehn@ymc.de

www.ymc.de