

Accessories for YMC Glass Columns



Accessories supplied with ECO and ECO^{PLUS} 10-15 mm ID

Connection Principle



Piston Thread: 1/4"-28G Ferrule Nut Tubing 1/16" (0.8 x 1.6)

Connection to
Column
1/4"-28 G



K.P316


Connection to
**Column/
System**
1/4"-28 G



JR-CFL-CB1KFS

JR-55050

Connection to
System
10-32 HPLC
(GE ÄKTA)



JR-58000

Connection to
System
M6
(GE ÄKTA)




K.P200N

K.P207

Accessories supplied with ECO and ECO^{PLUS} 25-50 mm ID

Connection Principle



Piston Thread: 1/4"-28G Ferrule Nut Tubing 1/8" (1.6 x 3.2)

Connection to
Column
1/4"-28 G



K.P316

Connection to
**Column/
System**
1/4"-28 G



JR-CFL-CB2KFS

JR-55051

Connection to
System
M6
(GE ÄKTA)



K.P300N

K.P307

Which Fitting Do I Have?

System

Fitting for the System
M6 (GE ÄKTA)



K.P207 1/16" Tubing
K.P307 1/8" Tubing

System

Fitting for the System
10-32 (GE ÄKTA)



JR-58000

Column

Lock Screw for the
Glass Column 1/4"-28 G



K.P316

Column/System

Fitting for the System/Glass Column
1/4"-28 G



JR-55050 1/16" Tubing
JR-55051 1/8" Tubing